



Covaris settings

Duty cycle	20%
Intensity	8
Cycles per burst	200
Cycle time	30 sec

Sample ID	ng/ μ l	260/280	260/230	Amount on gel (μ g)	Total amount (μ g/1.25 ml)
CTRL 3 min	191	1.83	1.52	3.83	239.30
LPS 3 min	203	1.84	1.73	4.06	253.49
CTRL 10 min	196	1.9	2.08	3.92	245.00
LPS 10 min	169	1.87	2.06	3.38	211.40