
IMPORTANT CUSTOMER INFORMATION

TruSeq™ Sample Prep Kit and Cluster Kit Compatibility

Dear Researcher,

Thank you for purchasing a TruSeq Sample Preparation Kit. Due to sequence changes in adapters and sequencing primers, libraries constructed with TruSeq Sample Preparation Kits might not work with older versions of Cluster Station and cBot cluster generation kits.

Sample Preparation Kit	TruSeq v2 on cBot TruSeq v5 on Cluster Station	v1 on cBot v4 on Cluster Station
TruSeq DNA/RNA/Exome	Compatible	Only compatible if index read primer and multiplexing read 2 primer are used
TruSeq Small RNA	Compatible	Not Compatible
Older DNA or mRNA	Compatible	Compatible
Small RNA v1.0/v1.5	Not Compatible	Compatible

For questions regarding product compatibility contact Illumina Technical Support.

Illumina, Inc.

9885 Towne Centre Drive
San Diego, CA 92121

Web: www.illumina.com
Email: techsupport@illumina.com
Phone: 800-809-ILMN (4566) (North America)
+1-858-202-ILMN (4566) (Outside North America)
Fax: 858-202-4766

